



VMRD, Inc.

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vtech@vmrd.com

Web site: <http://www.vmrd.com>

Certificate of Analysis

FA CONJUGATE DILUTING BUFFER

CATALOG NO.: 210-91-CB

VOLUME: 100 ml

LOT: V01109-052703

EXPIRATION: 27 May 2003

DESCRIPTION: Conjugate Diluting Buffer is PBS-0.09% azide with 1% BSA added.

INTENDED USE: Conjugate Diluting Buffer should be used for diluting concentrated conjugates for direct and indirect fluorescence.

QUALITY CONTROL METHOD: CDV, *Ehrlichia canis*, *Neospora caninum*, BVD, FIP-1, FIP-2, ERV/EHV, and CPV 10-well substrate slides tested with corresponding positive sera and anti-immunoglobulin FITC conjugates.

Specific Reaction: Dilute to same signal as previous lot.

Comments: Zero background.

STORAGE: Store at 4°C.

FOR RESEARCH AND INVESTIGATIONAL USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.