



**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)

Web site: <http://www.vmrd.com>

---

## Certificate of Analysis

### **FA MOUNTING FLUID**

**CATALOG NO.:** 210-92-MF

**VOLUME:** 10 ml

**LOT:** P060130-001

**EXPIRATION:** 8 February 2008

**DESCRIPTION:** FA Mounting Fluid pH 9.0-9.6 is made from a 50/50 mix of FA Rinse Buffer and glycerol.

**INTENDED USE:** FA Mounting Fluid is applied to FA substrate slides beneath the cover slip prior to viewing with fluorescence microscope.

**QUALITY CONTROL METHOD:** Compared to previous lot by direct FA using VMRD, Inc. IBR 10-well substrate slides (catalog no. SLD-IFA-IBR); IBR Direct FA Conjugate (catalog no. 210-69-IBR).

**Comments:** No background fluorescence.

**STORAGE:** FA Mounting Fluid should be stored at room temperature.

FOR *IN VITRO* LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.