

CERTIFICATE OF ANALYSIS

## Feline Infectious Peritonitis Virus (FIP)

FITC Direct Conjugate

Catalog No.:	CJ-F-FIP-10ML
Lot:	P141120-001
Expiration:	26 November 2018

**Quality Control Method:** Direct FA using FIP-1 2-well slide (catalog no. SLD-FAC-FIP1) and FIP-2 2-well slide (catalog no. SLD-FAC-FIP2).

**Specific Reaction:** 2-4+ fluorescence on the positive well and negative on the negative well, no background.

**Other Comments:** NA

**QC Manager:** *Torrey Beeson* Torrey Beeson

**QC Release Date:** 26 November 2014

**Pattern Of Fluorescence:** Multiple syncytia, plaques and some individual cells with granular cytoplasmic fluorescence.

**FA Cross-Reaction Screening:** The raw material has also been screened by indirect FA and has been found to react with *Toxoplasma gondii* at trace to 1+, but does not react with feline calicivirus (FCV), feline leukemia virus (FeLV), feline herpesvirus (FHV-1/FVR), feline panleukopenia virus (FPLV), reovirus (REO), rabies Recombinant Nucleoprotein (rNP), vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\COA (VMRD)\Conjugate\FIP\FIP - 10 ml\CJ-F-FIP-10ML\_P141120-001\_181126.docx

05 December 2014