

CERTIFICATE OF ANALYSIS

2A2

Monoclonal Antibody

<b>Catalog No. / Cell Line:</b>	2A2
<b>Isotype:</b>	IgG <sub>2a</sub>

**Specificity:**

Bovine Parainfluenza Virus Type 3, SF4 strain (69 kDa protein).

**Known Applications:**

Virus neutralization, hemagglutination inhibition.

**References:**

No specific reference at this time.

**Description and Handling:**

This monoclonal antibody is produced as mouse ascites fluid, clarified by centrifugation, and filtered through a 0.2 µM filter. The concentration is 1.0 mg/ml in phosphate-buffered saline, preserved with sodium azide.

**Storage:**

When the vial is stored at 2-8°C, it should be stable for one year.

**FOR IN VITRO LABORATORY USE ONLY.**