

**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815  
Fax: + 1 (509) 332-5356  
E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)  
Web site: <http://www.vmrd.com>

## Certificate of Analysis

### ***Toxoplasma gondii*** *Feline Negative Control Serum*

**CATALOG NO.:** 211-N-TOXO-FEL**VOLUME:** 1 ml**LOT:** P030227-003**EXPIRATION DATE:** 16 June 2008**DESCRIPTION:** Feline serum diluted in PBS, 1% BSA, 0.09% sodium azide.**INTENDED USE:** As negative control serum in detection of antibody to *Toxoplasma gondii* by indirect FA technique. This serum should be used undiluted to demonstrate non-specific fluorescence, if any.**QUALITY CONTROL METHOD:** Indirect FA using VMRD, Inc. *Toxoplasma* Substrate Slides (catalog no. SLD-IFA-TOXO); Anti-Feline IgG FITC Conjugate (catalog no. 055-10).**Specific Reaction:** Negative.**Other Reactions or Comments:** 2-3+ apical fluorescence.**PATTERN OF FLUORESCENCE:** None with our system.**STORAGE:** Store at 4-8°C.

FOR *IN VITRO* LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.