


Bartonella henselae IgM

FA Positive Control

Catalog No.:	PC-IFA-BH-M
Shelf Life:	4 years from date of qualification
Volume:	1 ml
Immunogen:	<i>Bartonella henselae</i>
Preservative:	0.09% Sodium Azide
Species of Origin:	Feline



Description:

Feline serum diluted in PBS, 10% DBS and 0.09% sodium azide as a preservative.

Intended Use:

As a positive control serum in detection of antibody to *Bartonella henselae* by indirect FA technique. This serum should be used undiluted to demonstrate positive fluorescence.

Storage:

Store at 2-7°C. DO NOT FREEZE! If control becomes cloudy it should be discarded.

References:

NA

Technical Data Sheet Version:

Version 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\TDS\Positive Controls\PC-IFA-BH-M TDS V1 CPM Active MKT.docx

25 September 2014

