

**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vmrd@vmrd.comWeb site: <http://www.vmrd.com>

Certificate of Analysis

Neospora caninum

Positive Control Serum

Canine Origin


CATALOG NO.: 211-P-NC-CAN

VOLUME: 1 ml

LOT: V01047-102707

EXPIRATION DATE: 27 October 2007

AGENT: *Neospora caninum*



DESCRIPTION: Canine sera diluted in PBS, 1% BSA, 0.09% sodium azide.

INTENDED USE: As control serum in detection of antibody to *Neospora caninum* by indirect FA technique. This control should be used undiluted to demonstrate positive fluorescence.

QUALITY CONTROL METHOD: Indirect FA using VMRD, Inc. *Neospora caninum* 12-well Substrate Slides (catalog no. 210-88-12-NC); Anti-Canine IgG FITC Conjugate (catalog no. 035-10); *Neospora caninum* Canine Negative Control Serum (catalog no. 211-N-NC-CAN).

Specific Reaction: 4+ positive fluorescence at neat.

Other Reactions or Comments: No background. IFA endpoint titer greater than 1:16.

PATTERN OF FLUORESCENCE: Peripheral fluorescence of *Neospora caninum* organisms. Apical (polar) staining is not considered positive. High titer test sera will fluoresce brighter than this positive control serum.

STORAGE: Store at 4-8°C. DO NOT FREEZE!

FOR IN VITRO LABORATORY USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.