

CERTIFICATE OF ANALYSIS

2D6

Monoclonal Antibody

Catalog No.:	2D6
Lot:	P181112-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using PRRSV-NA 12-well slide (catalog no. SLD-IFA-PRRS-NA), isotype control IgG2b and anti-murine IgG affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse Immunoglobulin IgG_{2b} kit.

Specific Reaction: 2 – 4+ signal at 10 µg/ml with an endpoint concentration less than 0.01 µg/ml. The isotype control was negative, no background. The RID test results showed a concentration of 1.194 mg/ml.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 05 March 2019

Pattern Of Fluorescence: Granular and diffuse fluorescence on cytoplasm, membrane, and occasional nuclear and cytoplasmic inclusions.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

11 March 2019