

CERTIFICATE OF ANALYSIS

2E10-H2

Monoclonal Antibody

Catalog No.:	2E10-H2
Lot:	P160817-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using CAV-1 in house slide, IgG1 isotype control, and anti-murine IgG affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse Immunoglobulins IgG1 kit.

Specific Reaction: 1-3+ signal at 10 µg/ml with an endpoint concentration greater than 1 µg/ml. The isotype control was negative with Tr-1+ background. The RID test results showed a concentration of 1.18 mg/ml.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 25 August 2016

Pattern of Fluorescence: Bright nuclear staining with moderate diffuse cytoplasmic fluorescence.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Monoclonals\2E10-H2\2E10-H2_P160817-001.docx

26 August 2016