

CERTIFICATE OF ANALYSIS

2E10-H2


Monoclonal Antibody

Catalog No.:	2E10-H2
Lot:	P230329-002
Expiration:	1 year when stored at 2-7 °C

Quality Control Method: Indirect FA using CAV-1 in house slide, IgG₁ isotype control, and anti-murine IgG affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-10ML). RID using in-house IgG₁ plate.

Specific Reaction: 2-3+ signal at 10 µg/ml with an endpoint concentration less than 0.1 µg/ml. The isotype control was negative. The RID test results showed a concentration of 1.19 mg/ml.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Bright nuclear staining with moderate diffuse cytoplasmic fluorescence.

Technical Data Sheet Version: 1