

CERTIFICATE OF ANALYSIS

348

Monoclonal Antibody

|              |                             |
|--------------|-----------------------------|
| Catalog No.: | 348                         |
| Lot:         | P180823-001                 |
| Expiration:  | 1 year when stored at 2-7°C |

**Quality Control Method:** Indirect FA using BVDV-1 (Modderman) 12-well (in house slide), isotype control IgG<sub>2b</sub>, and anti-murine affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse Immunoglobulin IgG<sub>2b</sub> kit.

**Specific Reaction:** 2 – 3+ signal at 10 µg/ml with no background, and an endpoint concentration less than 0.1 µg/ml. The isotype control was negative, no background. The RID test results showed a concentration of 0.987 mg/ml.

**Other Comments:** NA

**QC Manager:** Chris Cay



**QC Release Date:** 05 March 2019

**Pattern Of Fluorescence:** Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

**Technical Data Sheet Version:** Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

18 March 2019