

CERTIFICATE OF ANALYSIS

Canine Coombs Reagent

Coombs Reagent

Catalog No.:	392-5
Lot:	P160104-001
Expiration:	02 February 2018

Quality Control Method: Washed sheep red blood cells (SRBC) were sensitized with the canine Coombs positive control (catalog no. 372-2). The procedure is performed according to the Coombs positive control procedure. The Coombs reagent was tested according to the Coombs reagent procedure.

Specific Reaction: The Coombs reagent produced agglutination on sensitized SRBCs of 3+ agglutination at 1/2, 3+ at 1/4 and 2+ agglutination at 1/8. No agglutination was observed with sensitized SRBC in PBS without the Coombs reagent.

Other Comments: The subagglutinating dose for the canine positive control was determined to be 1/10. When using the positive control this dose may vary under your laboratory conditions and especially with your source of SRBCs. Therefore, we recommend you titer the positive control with your own SRBCs before using it. Please refer to Section A of the Coombs positive control procedure.

QC Manager:  **Will Thurman**

QC Release Date: 03 February 2016

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Coombs\Canine Coombs\Canine Coombs - 5 ml\392-5_P160104-001_180202.docx

10 February 2016