

CERTIFICATE OF ANALYSIS

3C9D11H11

Monoclonal Antibody

Catalog No.:	3C9D11H11
Lot:	P150408-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using PPV 12-well slide (catalog no. SLD-IFA-PPV), isotype control IgG₁ (1.0mg/ml), and anti-murine FITC affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse RID kit IgG₁.

Specific Reaction: T-2+ fluorescence at 10 µg/ml with an endpoint concentration less than 1.0 µg/ml. The isotype control was negative with no background. The RID test results showed a concentration of 1.0 mg/ml.

Other Comments: NA

QC Manager: *Torrey Beeson* Torrey Beeson

QC Release Date: 20 April 2015

Pattern Of Fluorescence: Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\COA (VMRD)\Monoclonals\3C9D11H11\3C9D11H11_P150408-001.docx

22 April 2015