

**CERTIFICATE OF ANALYSIS**

3C9D11H11

Monoclonal Antibody

Catalog No.:	3C9D11H11
Lot:	P170424-001
Expiration:	1 year when stored at 2-7°C

**Quality Control Method:** Indirect FA using PPV 12-well slide (catalog no. SLD-IFA-PPV), isotype control IgG1 (1.0mg/ml), and anti-murine FITC affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse RID kit IgG1.

**Specific Reaction:** 2-4+ signal at 10 µg/ml with no background. Endpoint titer was found to be less than 1.0 µg/ml. RID testing showed a concentration of 0.831 mg/ml.

**Other Comments:** NA

**QC Manager:**  Will Thurman

**QC Release Date:** 10 May 2017

**Pattern of Fluorescence:** Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Monoclonal\3C9D11H11\3C9D11H11\_P170424-001.docx

02 June 2017