

**CERTIFICATE OF ANALYSIS**

## 3G9F3 Pseudorabies Virus (PRV) gIII

Monoclonal Antibody

Catalog No.:	3G9F3
Lot:	P171030-001
Expiration:	1 year when stored at 2-7°C

**Quality Control Method:** Indirect FA using PRV 12-well slide (in house), isotype control IgG<sub>2b</sub> and anti-murine affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Use mouse Immunoglobulins kit IgG<sub>2b</sub>.

**Specific Reaction:** 1 – 3+ fluorescence at 10 µg/ml with an endpoint concentration less than 0.01 µg/ml. The isotype control was negative with no background. The RID test results showed a concentration of 1.091 mg/ml.

**Other Comments:** NA

**QC Manager:**  Will Thurman

**QC Release Date:** 14 November 2017

**Pattern Of Fluorescence:** Plaques of infected cells exhibiting granular cytoplasmic fluorescence with occasional nuclear staining.

**Technical Data Sheet Version:** Version #2

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Master Documents\COA (VMRD)\Monoclonal\3G9F3\3G9F3\_P171030-001.docx

13 December 2017