

CERTIFICATE OF ANALYSIS

## 3G9F3 Pseudorabies Virus (PRV)

Monoclonal Antibody

Catalog No.:	3G9F3
Lot:	P210621-004
Expiration:	1 year when stored at 2-7 °C

**Quality Control Method:** Indirect FA using PRV 12-well slide (in house), isotype control IgG2b and anti-murine affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Isotype concentration determined by RID.

**Specific Reaction:** 2-4+ signal at 10 µg/ml with an endpoint concentration less than 0.1 µg/ml. The isotype control was negative with no background. The RID test results showed a concentration of 0.889 mg/ml.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Plaques of infected cells exhibiting granular cytoplasmic fluorescence with occasional nuclear staining.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

08 October 2021