

CERTIFICATE OF ANALYSIS

3.12F1

Monoclonal Antibody

| | |
|--------------|-----------------------------|
| Catalog No.: | 3.12F1 |
| Lot: | P180820-008 |
| Expiration: | 1 year when stored at 2-7°C |

Quality Control Method: Indirect FA using BVDV 12-well slide (catalog no. SLD-IFA-BVD), isotype control IgG₁ and anti-murine affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse RID kit IgG₁.

Specific Reaction: 1 – 3+ signal with no background at 10 µg/ml with an endpoint concentration less than 10 µg/ml. The isotype control was negative, no background. The RID test results showed a concentration of 1.075 mg/ml.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 03 October 2018

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

04 October 2018