

**CERTIFICATE OF ANALYSIS**

3.12F1

Monoclonal Antibody

Catalog No.:	3.12F1
Lot:	P221007-001
Expiration:	1 year when stored at 2-7 °C

**Quality Control Method:** Indirect FA using BVDV in-house 12-well slide, isotype control IgG<sub>1</sub> and anti-murine affinity purified conjugate (catalog no. CJ-F-MURG-AP-10ML). RID with in-house kit.

**Specific Reaction:** 1–3+ signal with no background at 10 µg/ml with an endpoint concentration less than 10 µg/ml. The isotype control was negative with trace–1+ background. The RID test results showed a concentration of 0.921 mg/ml.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

18 November 2022