

CERTIFICATE OF ANALYSIS

Feline Coombs Reagent

Coombs Reagent

Catalog No.:	592-2
Lot:	P150820-001
Expiration:	22 October 2017

Quality Control Method: Antibody to feline IgG and IgM is evaluated using immunodiffusion gel. Specificity is evaluated by testing the reagent according to the Coombs Reagent Procedure against erythrocytes from a number of clinically normal felines.

Other Comments: N/A

QC Manager:  **Will Thurman**

QC Release Date: 22 October 2015

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COAs and TDS\COA (VMRD)\Coombs\Feline Coombs\592-2_P150820-001_171022.docx

11 November 2015