

CERTIFICATE OF ANALYSIS

Feline Coombs Reagent

Coombs Reagent

Catalog No.:	592-2
Lot:	P170301-002
Expiration:	21 June 2019

Quality Control Method: Antibody to feline IgG and IgM is evaluated using immunodiffusion gel. Specificity is evaluated by testing the reagent according to the Coombs Reagent Procedure against erythrocytes from a number of clinically normal felines.

Other Comments: NA

QC Manager:  **Will Thurman**

QC Release Date: 21 June 2017

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Coombs\Feline Coombs\592-2_P170301-002_190622.docx

23 June 2017