

CERTIFICATE OF ANALYSIS

Feline Coombs Reagent

Coombs Reagent

Catalog No.:	592-2
Lot:	P180214-001
Expiration:	04 May 2020

Quality Control Method: Antibody to feline IgG and IgM is evaluated using immunodiffusion gel. Specificity is evaluated by testing the reagent according to the Coombs Reagent Procedure against erythrocytes from a number of clinically normal felines.

Other Comments: NA

Production Manager:



Erin Garrison

QC Release Date: 04 May 2018

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

04 May 2018