

CERTIFICATE OF ANALYSIS

Feline Coombs Reagent


Coombs Reagent

Catalog No.:	592-2
Lot:	P210504-001
Expiration:	19 December 2025

Quality Control Method: Antibody to feline IgG and IgM is evaluated using immunodiffusion gel. Specificity is evaluated by testing the reagent according to the Coombs Reagent Procedure against erythrocytes from a number of clinically normal felines.

Other Comments: NA

Production Manager:



Amanda Bronowski

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

22 December 2023