

CERTIFICATE OF ANALYSIS

## Feline Coombs Reagent


Coombs Reagent

Catalog No.:	592-2
Lot:	P220212-001
Expiration:	19 December 2025

**Quality Control Method:** Antibody to feline IgG and IgM is evaluated using immunodiffusion gel. Specificity is evaluated by testing the reagent according to the Coombs Reagent Procedure against erythrocytes from a number of clinically normal felines.

**Other Comments:** NA

**Production Manager:**



Amanda Bronowski

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

22 December 2023