

CERTIFICATE OF ANALYSIS

Bovine Leukemia Virus (BLV) MAb Ascites gp51 D-D'
IgG¹ Isotype

Monoclonal Antibody

Catalog No.:	BLV2
Lot:	P171212-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using BLV 12-well slide (catalog no. SLD-IFA-BLV), isotype control IgG¹, and anti-murine IgG affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Use mouse immunoglobulins kit IgG¹.

Specific Reaction: Tr-2+ signal at 10 µg/ml with an endpoint concentration less than <10 µg/ml. The isotype control was negative, with trace background. The RID test results showed a concentration of 0.884 mg/ml.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 25 January 2018

Pattern Of Fluorescence: Smooth, undifferentiated cytoplasmic fluorescence, cytoplasmic punctate inclusions in single cells and small syncytia (3-6 nuclei).

Technical Data Sheet Version: Version # 2

FOR IN VITRO LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Monoclonal\BLV2\BLV2_P171212-001_190125.docx

12 February 2018