

CERTIFICATE OF ANALYSIS

Canine Adenovirus (CAV)

FITC Conjugate

Catalog No.:	CJ-F-CAV-10ML
Lot:	P161122-001
Expiration:	01 December 2020

Quality Control Method: Direct FA on CAV-2 2 well slide (SLD-FAC-CAV).

Specific Reaction: 2-3+ signal with Tr-1+ background.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 02 December 2016

Pattern of Fluorescence: Mostly large, single nuclear inclusion and some cytoplasmic fluorescence in a few infected cells.

FA Cross-Reaction Screening: The raw material used in the production of this conjugate has been screened by FA and has been found to react with Canine Adenovirus Type1 (CAV-1) at a 1/200 dilution with a Tr-2+ signal. This is equivalent to this product at neat. The raw material does not react with Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Conjugate\CAV\CAV - 10 ml\CJ-F-CAV-10ML_P161122-001_201201.docx

07 December 2016