

CERTIFICATE OF ANALYSIS

Porcine Adenovirus (PAV)

FITC Conjugate

Catalog No.:	CJ-F-PAV-10ML
Lot:	P150904-001
Expiration:	16 October 2019

Quality Control Method: Must produce at least 3+ specific signal and equal to or less than 1+ background.

Specific Reaction: 2-4+ signal with trace background on the positive well and negative with trace background on the negative well.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 16 October 2015

Pattern Of Fluorescence: Granular nuclear membrane fluorescence on individual cells or plaques of cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Hemagglutinating Encephalomyelitis Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Transmissible Gastroenteritis Virus (TGEV), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COAs and TDS\COA (VMRD)\Conjugate\PAV\PAV - 10 ml\CJ-F-PAV-10ML_P150904-001_191016.docx

19 October 2015