

CERTIFICATE OF ANALYSIS

Porcine Adenovirus (PAV)


FITC Conjugate

Catalog No.:	CJ-F-PAV-10ML
Lot:	P220527-001
Expiration:	19 September 2030

Quality Control Method: Direct FA on PAV 2-well control slides (in house).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well. No background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Granular nuclear membrane fluorescence on individual cells or plaques of cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Hemagglutinating Encephalomyelitis Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Transmissible Gastroenteritis Virus (TGEV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

23 September 2022