

CERTIFICATE OF ANALYSIS

# Porcine Hemagglutinating Encephalomyelitis Virus (PHEV)

FITC Conjugate

Catalog No.:	CJ-F-PHEV-10ML
Lot:	P190805-001
Expiration:	01 August 2031

Quality Control Method: Direct FA using PHEV 12-well slide.

Specific Reaction: 3-4+ signal with no background.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic and membrane fluorescence on individual cells and some syncytia.

**FA Cross-Reaction Screening:** The raw material has also been screened by FA and does not react with Porcine Adenovirus (PAV), Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Transmissible Gastroenteritis Virus (TGEV), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

03 August 2023