

CERTIFICATE OF ANALYSIS

Porcine Hemagglutinating Encephalomyelitis Virus (PHEV)

FITC Conjugate

Catalog No.:	CJ-F-PHEV-10ML
Lot:	P230710-001
Expiration:	19 July 2031

Quality Control Method: Direct FA using a PHEV in-house 12-well slide

Specific Reaction: 3–4+ signal with no background.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic and membrane fluorescence on individual cells and some syncytia.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Porcine Adenovirus (PAV), Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Transmissible Gastroenteritis Virus (TGEV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

24 July 2023