

**CERTIFICATE OF ANALYSIS**

# Porcine Parvovirus (PPV)

FITC Conjugate

Catalog No.:	CJ-F-PPV-10ML
Lot:	P141210-002
Expiration:	29 October 2022

Quality Control Method: Direct FA using PPV 2-well slide (catalog no. SLD-FAC-PPV).

Specific Reaction: 1-3+ signal on the positive well, no background. Negative on the negative well, no background.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 30 October 2018

Pattern Of Fluorescence: Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

**FA Cross-Reaction Screening:** The raw material has also been screened by indirect FA and has been found to react with porcine cytomegalovirus (PCMV) 3+ at 1/50, but does not react with porcine adenovirus (PAV), porcine circovirus type 1 and 2 (PCV-1 and 2), porcine hemagglutinating encephalomyelitis virus (PHEV), porcine reproductive and respiratory syndrome virus North American and European strains (PRRSV), transmissible gastroenteritis virus (TGEV), reovirus (REO), rabies Recombinant Nucleoprotein (rNP), vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

5 November 2018