

CERTIFICATE OF ANALYSIS

Transmissible Gastroenteritis Virus (TGEV)

FITC Conjugate

Catalog No.:	CJ-F-TGE-10ML
Lot:	P160713-002
Expiration:	12 October 2020

Quality Control Method: Direct FA using TGEV 2-well slides (catalog no. SLD-FAC-TGE).

Specific Reaction: 2-4+ signal on the positive well, negative on the negative well with trace background.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 13 October 2016

Pattern of Fluorescence: Moderate to large syncytia with granular and smooth cytoplasmic fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Porcine Adenovirus (PAV), Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Hemagglutinating Encephalomyelitis Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Conjugate\TGE\TGE - 10 ml\CJ-F-TGE-10ML_P160713-002_201012.docx

19 October 2016