

CERTIFICATE OF ANALYSIS

Transmissible Gastroenteritis Virus (TGEV)


FITC Conjugate

Catalog No.:	CJ-F-TGE-10ML
Lot:	P210511-001
Expiration:	14 May 2025

Quality Control Method: Direct FA using TGEV 2-well slides (catalog no. SLD-FAC-TGE).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well, trace background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Moderate to large syncytia with granular and smooth cytoplasmic fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Porcine Adenovirus (PAV), Porcine Circovirus Type 1 and 2 (PCV-1 and 2), Porcine Hemagglutinating Encephalomyelitis Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European strains, Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

19 May 2021