

CERTIFICATE OF ANALYSIS

CPI-A-CA

Monoclonal Antibody

Catalog No.:	CPI-A-CA
Lot:	P150318-002
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using CPI-2 12-well slide (catalog no. SLD-IFA-CPI), isotype control IgG₁ and anti-murine IgG FITC affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse RID kit IgG₁.

Specific Reaction: 3-4+ fluorescence at 10 µg/ml with an endpoint concentration less than 1.0 µg/ml. The isotype control was negative with no background. The RID test results showed a concentration of 1.04 mg/ml.

Other Comments: NA

QC Manager: *Torrey Beeson* Torrey Beeson

QC Release Date: 24 March 2015

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.