

**CERTIFICATE OF ANALYSIS**

## CPI-A-CA

Monoclonal Antibody

Catalog No.:	CPI-A-CA
Lot:	P151015-001
Expiration:	1 year when stored at 2-7°C

**Quality Control Method:** Indirect FA using CPI-2 12-well slide (catalog no. SLD-IFA-CPI), isotype control IgG1 and anti-murine IgG FITC affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse RID kit IgG<sub>1</sub>.

**Specific Reaction:** 1-3+ signal at 10 µg/ml, endpoint reached at 1 µg/ml. Isotype control was negative. RID concentration was 0.84 mg/ml.

**Other Comments:** NA

**QC Manager:**  Will Thurman

**QC Release Date:** 20 October 2015

**Pattern Of Fluorescence:** Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

Technical Data Sheet Version: Version #1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality\COAs and TDS\COA (VMRD)\Monoclonals\CPI-A-CA\CPI-A-CA\_P151015-001.docx

23 October 2015