

CERTIFICATE OF ANALYSIS

CPI-A-CA Canine Parainfluenza Virus Type 2 (CPI-2)

Monoclonal Antibody

Catalog No.:	CPI-A-CA
Lot:	P160406-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using CPI-2 12-well slide (catalog no. SLD-IFA-CPI), isotype control IgG1 and anti-murine IgG FITC affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse IgG1 kit.

Specific Reaction: 2-3+ signal at 10 µg/ml dilution. Endpoint reached at 1 µg/ml. The isotype control was negative with no background. The RID test results showed a concentration of 1.03 mg/ml.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 12 April 2016

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Monoclonals\CPI-A-CA\CPI-A-CA_P160406-001.docx

14 April 2016