

CERTIFICATE OF ANALYSIS

D89

Monoclonal Antibody

Catalog No.:	D89
Lot:	P180430-003
Expiration:	1 year when stored at 2-7°C

Quality Control Method: Indirect FA using BVDV-1 (Modderman) 12-well (in house slide), isotype control IgG_{2a} and anti-murine affinity purified conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). Mouse Immunoglobulin IgG_{2a} kit.

Specific Reaction: 2 – 3+ signal at 10 µg/ml with an endpoint concentration less than 10 µg/ml. The isotype control was negative, no background. The RID test results showed a concentration of 0.985 mg/ml.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 21 June 2018

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

21 June 2018