

CERTIFICATE OF ANALYSIS

D89

Monoclonal Antibody

Catalog No.:	D89
Lot:	P211103-001
Expiration:	1 year when stored at 2-7 °C

Quality Control Method: Indirect FA using BVDV-1 (Modderman) 12-well (in house slide), isotype control IgG_{2a} and anti-murine affinity purified conjugate (catalog no. CJ-F-MURG-AP-10ML). RID performed with in-house RID kit.

Specific Reaction: 2-4+ signal at 10 µg/ml with an endpoint concentration less than 1.0 µg/ml. The isotype control was negative with trace background. The RID test results showed a concentration of 1.006 mg/ml.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: 1