

CERTIFICATE OF ANALYSIS

F99/97.6.1


Monoclonal Antibody

Catalog No.:	F99/97.6.1
Lot:	P210428-001
Expiration:	1 year when stored at 2-7 °C

Quality Control Method: The concentration was tested in RID with a mouse Immunoglobulins kit IgG1.

Specific Reaction: The RID test results showed a concentration of 1.049 mg/ml.

Other Comments: NA

QC Manager: 

David Hyatt

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

05 November 2021