

CERTIFICATE OF ANALYSIS

F99/97.6.1

Monoclonal Antibody

Catalog No.:	F99/97.6.1
Lot:	P190628-001
Expiration:	1 year when stored at 2-7°C

Quality Control Method: F99/97.6.1 cell max was evaluated by immunohistochemistry (IHC) of strong and weak tissues from scrapie-infected sheep, and a negative tissue sample from a sheep with no known exposure to scrapie. The antibody was diluted to 1.0 µg/ml and run according to the kit insert for VMRD Bovine Spongiform Encephalopathy Antigen Test Kit, Immunohistochemistry (catalog no. 298). The concentration was tested in RID with a mouse Immunoglobulins kit IgG₁.

Specific Reaction: There was no staining on the tissue from the sheep with no exposure to scrapie. Staining for the scrapie-infected sheep was 1+ on the strong positive and 1-2+ on the weak positive. RID results showed a concentration of 1.145 mg/ml.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 20 May 2019

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

21 May 2019