

CERTIFICATE OF ANALYSIS

FA Conjugate Diluting Buffer

Catalog No.:	FACDB-100ML
Lot:	P180612-002
Expiration:	22 June 2022

Quality Control Method: Direct FA using 1/50 dilution of 50X Conjugates and the corresponding slides on each of the following VMRD systems: Bovine PI-3 (catalog no. CJ-F-PI3-50X, and SLD-IFA-PI3), Canine CDV (catalog no. CJ-F-CDV-50X, and SLD-IFA-CDV), Equine/Llama ERV (catalog no. CJ-F-ERV-50X, and SLD-IFA-ERV), Feline FVR (catalog no. CJ-F-FVR-50X, and SLD-IFA-FVR), and Porcine TGE (catalog no. CJ-F-TGE-50X, and SLD-IFA-TGE).

Specific Reaction: 1-4+ signal for all systems at 1/50 with no additional background compared to the previous lot. CDB alone produced no background and no signal with all systems.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 25 June 2018

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

26 June 2018