

**CERTIFICATE OF ANALYSIS**

## FA Mounting Fluid

Catalog No.:	FAMF-10ML
Lot:	P150602-001
Expiration:	08 June 2019

**Quality Control Method:** Direct FA using BHV-1/IBR 2-well slides (catalog no. SLD-FAC-IBR) and FITC conjugate (catalog no. CJ-F-IBR-1ML or 10ML).

**Specific Reaction:** There is no difference between the reference and the QC lot of mounting fluid.

pH: 9.24

**Other Comments:** NA

**QC Manager:** *Torrey Beeson* Torrey Beeson

**QC Release Date:** 10 June 2015

**Pattern Of Fluorescence:** The pattern of fluorescence will vary depending on what system was used.

**Technical Data Sheet Version:** 1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality VMRD\COAs and TDS\COA (VMRD)\Buffers\Mounting Fluid\FAMF-10ML\_P150602-001\_190608.docx

12 June 2015