

CERTIFICATE OF ANALYSIS

## FA Mounting Fluid

Catalog No.:	FAMF-10ML
Lot:	P170517-003
Expiration:	24 May 2021

Quality Control Method: Direct FA using BHV-1/IBR 2-well slides (catalog no. SLD-FAC-IBR) and FITC conjugate (catalog number: CJ-F-IBR-1ML or 10ML).

Specific Reaction: No difference in signal and background noted between the new lot and the previous lot.

pH: 9.54

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 25 May 2017

Pattern of Fluorescence: The pattern of fluorescence will vary depending on what system was used.

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Buffers\Mounting Fluid\FAMF-10ML\_P170517-003\_210524.docx

02 June 2017