

CERTIFICATE OF ANALYSIS

4X FA Powdered Rinse Buffer


Catalog No.:	FARB-4X
Lot:	P220130-001
Expiration:	10 June 2026

Quality Control Method: Reconstituted in 1 liter of deionized water to make 4X concentrate. Confirmed pH of reconstituted solution to be between 9.0 and 9.5 with pH meter. Diluted 4X solution to 1X and used for direct FA with FIP conjugate (catalog no. CJ-F-FIP-10ML) and FIP2 2-well control slide (catalog no. SLD-FAC-FIP2).

Specific Reaction: Reconstituted 4X and 1X rinse buffer was clear. QC rinse buffer produced results similar to reference rinse buffer.

pH: 9.09

Other Comments: NA

QC Manager: 

David Hyatt

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

15 June 2022