

**CERTIFICATE OF ANALYSIS**

# Canine Parainfluenza Virus Type 2 (CPI-2)

Antiserum

Catalog No.:	PAB-CPI
Lot:	P150309-002
Expiration:	20 March 2019

**Quality Control Method:** Indirect FA using CPI-2 12-well slide (catalog no. SLD-IFA-CPI), and anti-porcine IgG FITC affinity purified conjugate (catalog no. CJ-F-PORG-AP-1ML or 10ML).

**Specific Reaction:** 1-4+ fluorescence at 1/400 with an endpoint titer equal to or greater than 1/10,000.

**Other Comments:** NA

**QC Manager:** *Torrey Beeson* Torrey Beeson

**QC Release Date:** 20 March 2015

**Pattern Of Fluorescence:** Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

**FA Cross-Reaction Screening:** The Antiserum has also been screened by Indirect FA and does not react with Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parvovirus (CPV), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality VMRD\COAs and TDS\COA (VMRD)\Antisera\CPI\PAB-CPI\_P150309-002\_190320.docx

26 March 2015