

CERTIFICATE OF ANALYSIS

Canine Parainfluenza Virus Type 2 (CPI-2)

Antiserum

Catalog No.:	PAB-CPI
Lot:	P181101-001
Expiration:	22 February 2023

Quality Control Method: Indirect FA using CPI-2 12-well slide (catalog no. SLD-IFA-CPI), and anti-porcine IgG FITC affinity purified conjugate (catalog no. CJ-F-PORG-AP-1ML or 10ML).

Specific Reaction: 2 - 4+ signal at 1/400 with an endpoint titer equal to or greater than 1/10,000.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 22 February 2019

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

FA Cross-Reaction Screening: The Antiserum has also been screened by Indirect FA and does not react with Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parvovirus (CPV), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

12 March 2019