

CERTIFICATE OF ANALYSIS

Canine Parainfluenza Virus Type 2 (CPI-2)

Antiserum

| | |
|--------------|------------------|
| Catalog No.: | PAB-CPI |
| Lot: | P210216-001 |
| Expiration: | 22 February 2025 |

Quality Control Method: Indirect FA using CPI-2 2-well slide (catalog no. SLD-FAC-CPI), and anti-porcine IgG FITC affinity purified conjugate (catalog no. CJ-F-PORG-AP-10ML)

Specific Reaction: 2–4+ signal with the antisera at 1/400 dilution, no background. Endpoint titer is 1/10,000.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells

FA Cross-Reaction Screening: The Antiserum has also been screened by Indirect FA and does not react with Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parvovirus (CPV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

05 March 2021