

CERTIFICATE OF ANALYSIS

Canine Parvovirus (CPV)


IgG FA Positive Control

Catalog No.:	PC-IFA-CPV-G
Lot:	P230611-001
Expiration:	22 June 2031

Quality Control Method: Indirect FA using CPV 12-well slide (catalog no. SLD-IFA-CPV) and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 3–4+ signal with the control at neat. No background. Endpoint signal was reached at 1/16 dilution.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Single cells with cytoplasmic and nuclear fluorescence.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

26 June 2023