


CERTIFICATE OF ANALYSIS

Equine Infectious Anemia Virus (EIAV), ELISA Medium Positive

Reference Serum

Catalog No.:	RS-EIA-EM-0.5ML
Size:	0.5 ml
Lot:	P120601-002
Expiration:	01 February 2025
Species of Origin:	Equine



Description:

Equine serum containing a level of antibody to Equine Infectious Anemia Virus that gives a medium positive reaction in the VMRD ELISA and a weak positive in the VMRD AGID.

Sample	OD Values
QC serum	0.341
+ Control	0.204

Quality Control Method:

Sample was tested in the VMRD EIA ELISA v2 (catalog no. 5515.01-1) and AGID (catalog no. 400-200) according to each kit insert.

ELISA Results: 0.341 OD

AGID Results: Weak positive

Other Information:

The sample's expected OD range in the VMRD EIA ELISA is > 0.25 and < 0.75. It is also usually weakly positive by AGID.

Intended Use:

This product is intended as a reference sample for quality assurance of EIA ELISA tests.

Storage:

This reference serum is provided in liquid form and should be stored at 2-7 °C. If serum becomes cloudy, it should be discarded. This serum contains 10 ppm ProClin 300 as a preservative.

References: NA

FOR *IN VITRO* LABORATORY USE ONLY.

03 February 2023