


CERTIFICATE OF ANALYSIS

EQUINE INFECTIOUS ANEMIA, ELISA STRONG POSITIVE

Reference Serum

Catalog No.:	RS-EIA-ES-0.5ML
Size:	0.5 ml
Lot:	P130322-002
Expiration:	28 April 2017
Species of Origin:	Equine



Description:

Equine serum containing a level of antibody to Equine Infectious Anemia Virus that gives a strong positive reaction in the VMRD ELISA and a positive in the VMRD AGID.

Sample	OD Values		Average OD
QC serum	0.664	0.712	0.688
+ Control	0.166	0.159	0.163
	0.165	0.162	

Quality Control Method:

Sample was tested in the VMRD EIA ELISA (catalog no. 290-1) and AGID (catalog no. 400-200) according to each kit insert.

ELISA Results: Average OD is 0.688

AGID Results: Positive

Other Information:

The sample's expected OD range in the VMRD EIA ELISA is > 0.5. It is also usually clearly positive by AGID.

Intended Use:

This product is intended as a reference sample for quality assurance of EIA ELISA tests.

Storage:

This reference serum is provided in liquid form and should be stored at 2-7°C. If serum becomes cloudy, it should be discarded. This serum contains 10 ppm ProClin 300 as a preservative.

References: NA

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\COA (VMRD)\Reference Samples\RS-EIA-ES-0.5ML_P130322-002_170428.docx
30 April 2015