


CERTIFICATE OF ANALYSIS

Equine Infectious Anemia Virus (EIAV), ELISA Weak Positive

Reference Serum

Catalog No.:	RS-EIA-EW-0.5ML
Size:	0.5 ml
Lot:	P120601-001
Expiration:	23 March 2023
Species of Origin:	Equine



Description:

Equine serum containing a level of antibody to Equine Infectious Anemia Virus that gives a weak positive reaction in the VMRD ELISA and a negative read in the VMRD AGID.

Sample	OD Values
QC serum	0.251
+ Control	0.221

Quality Control Method:

Sample was tested in the VMRD EIA ELISA v2 (catalog no. 5515.01-1) and AGID (catalog no. 400-200) according to each kit insert.

ELISA Results: OD is 0.251

AGID Results: Negative

Other Information:

The sample's expected OD range in the VMRD EIA ELISA is > 0.1 and < 0.4. It is also usually negative by AGID.

Intended Use:

This product is intended as a reference sample for quality assurance of EIA ELISA tests.

Storage:

This reference serum is provided in liquid form and should be stored at 2-7 °C. If serum becomes cloudy, it should be discarded. This serum contains 10 ppm ProClin 300 as a preservative.

References: NA

FOR IN VITRO LABORATORY USE ONLY.

25 March 2021